PERPOS™ PLS System Surgical Technique

This brochure is intended to demonstrate general surgical technique. Interventional Spine, Inc. as the manufacturer does not recommend this or any other surgical technique for use on a specific patient. The surgeon who performs the implant procedure is responsible for determining and utilizing the appropriate technique for implanting the device in each individual patient. Interventional Spine, Inc. is not responsible for selection of the appropriate surgical technique to be utilized for an individual patient.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. Interventional Spine®, Inc. is certified to IS EN ISO 13485:2003. The products have been assessed in conjunction with the Notified Body as applicable, and are considered to meet the Essential Requirements and so bear the CE Marking of Conformity.

As of date of print, Interventional Spine has several issued and pending U.S. patents. Interventional Spine®, PERPOS™ PLS System, BONE-LOK®, CLASP®, and Teleport® are all marks registered with the U.S. Patent and Trademark Office.

Interventional Spine, Inc.
13700 Alton Parkway, Suite 160
Irvine, CA 92618
Ph. (949) 472-0006
Fax (949) 472-0016
(800) 497-0484
www.i-spineinc.com

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

Step 2

Entry Point

The entry point is the location on the patient's skin where the Access Needle (9057) is inserted. Typically the entry point will be at the spinous process of the vertebral body one or two levels above the level to be treated (i.e. L3 for a L5/S1). This is a general guide, but patient anatomy will dictate actual position.

Step 3

Trajectory

Using Lateral and AP fluoroscopy images and the anatomy, the trajectory of the implant is determined, as shown in Fig 3 and Fig 4 below. The AP trajectory is medial to lateral.

Step 3 (cont.)

Trajectory

One technique to help see the proper trajectory is to place a K-wire near or on the patient while the fluoroscopy images are being taken. Adjusting the position of the wire will determine the best trajectory. This information is used to determine the Entry Point on the patient's skin and Start Point on the facet.

A small mid-line stab incision is made to accommodate the Access Needle. The Access Needle is inserted, with the bevel facing up (to prevent the Access Needle from "slipping" off of the facet), through the skin, along the proper trajectory until the tip is at the Start Point on the bony anatomy. A series of AP and Lateral views may be needed to confirm precise positioning. At this point the bevel of the access needle can be rotated to guide the K-wire. The wire will favor the long side of the beveled tip. This can be used in the AP or Lateral views to improve the trajectory.


Figure 3 - AP View Trajectory


Figure 4 - Lateral View Trajectory


Figure 5 - Access Needle

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Step 4 PERPOS™ K-Wire

With the inner stylet removed, the .062 K-Wire (6072-00) is passed through the Access Needle and driven across the facet joint, through the center of the pedicle, to (approximately) mid-portion of the inferior vertebral body (Fig. 6). This should be accompanied by a series of fluoroscopy images to insure the K-wire tip reaches the appropriate position and that the trajectory does not change as the wire goes across the facet joint into the center of the pedicle. The Access Needle is then removed while maintaining wire placement.


Figure 6 - K-wire

Step 5 Dilation

Once the K-Wire has been placed, make a 15mm incision through skin and fascia. The fascia needs to be incised inferiorly to the skin incision to accommodate Interventions Spine, Inc.'s 11mm ID Teleport® (6086). (Fig 7 & 8).

The Teleport® is placed over the K-wire and advanced until the first retractor comes in contact with the facet. This should be confirmed in the lateral view.


Figure 7 & 8 - Dilating with the Interventions Spine, Inc. 11mm ID Teleport® (6086)

Step 5 (cont.) Dilation

Pull the first button (labeled “Pull 1”) of the Teleport®, and then apply forward pressure to advance the second tissue retractor down to the facet. Pull the second button of the Teleport® (labeled “Pull 2”) and then apply forward pressure to advance the third retractor down to the facet. Use a lateral fluoroscopic image to verify proper placement of the Teleport®. After confirmation of the final placement of the outer dilator of the Teleport®, hold the flat distal handle in place, remove the proximal handle. The first and second dilator will come out and remain as one unit.

For procedures where the Teleport® is utilized more than once, remove the instrument from the initial surgical location. Reassemble the Teleport® by pulling the first button (labeled “Pull 1”), releasing the first retractor from the second retractor until “locked” in place. Place the third retractor over the first and second retractors, pull the first button (labeled “Pull 2”), and “lock into place. The Teleport® is now ready for use. Repeat the instructions above for the placement of the next BONE-LOK® PLS Implant (LSW-45-3040).*

*Note: The washer on the 4.5mm BONE-LOK® PLS Implant is 10.5mm in diameter - The selection of a proper retraction system (Minimum 11mm ID) is necessary to be able to accommodate the device.

Step 6 PERPOS™ 2-in-1 Drill

The cannulated 2-in-1 Drill (6070) is passed over the K-Wire. The drill should be connected to a cannulated power drill with an AO style quick connect or Jacob’s chuck (Fig. 9). The drill is advanced under power until the depth limiting feature (30mm “positive stop”) contacts the bone. This is confirmed through fluoroscopy.

Figure 9 - 2 in 1 Drill

The fluoroscopy image is saved and “flipped” to the second fluoroscopy screen. The tip of the drill is 30mm from the sphere - This image can be used as a reference for gauging the depth of the tap and placement of the BONE-LOK® PLS Implant.

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Step 7

PERPOS™ Tap

The cannulated Tap (6111) is driven by hand (Fig 10) into the bone until the tip reaches the appropriate depth (slightly past the distal tip of the drill, identified on the saved image from the second monitor), which is determined and verified by fluoroscopy. Once the appropriate depth is achieved, the tap is removed.

Figure 10 - Tap

Step 8

Implant Placement

The 4.5mm BONE-LOK® PLS Implant is fully seated onto the Driver (6112) and care is taken not to remove the plastic clip until implant is placed over the K-Wire. The anti-compression clip is then removed (Fig 11) and the device is advanced to, and then into, the bone. The placement of the implant should be with a continuous forward clockwise motion, using a series of fluoroscopy images to insure proper depth. Do not “pull back” on the Driver when inserting the implant – This may disengage the implant from the Driver, which could cause the premature advancement of the compression collar of the implant.

Figure 11 - BONE-LOK® PLS Implant

Step 8 (cont.) Implant Placement

The 4.5mm BONE-LOK® PLS Implant is inserted until the distal tip is at the desired location (as determined by the surgeon), leaving the proximal head of the implant slightly “proud” to allow for compression. After the proper position of the distal tip has been established, the driver is removed by pulling straight back.

(Note: Step 9 is used to fully seat the implant to the bone and gain compression of the facet joint.)

Step 9 Compression of Implant

The Compression Tool (6113-00) is passed over the K-Wire and over the Pull Pin (Fig 12).


Figure 12 - Compression

The device is designed to accommodate a range of 30 to 40 mm. This leaves a maximum of 10mm of available compression, but the actual amount of compression needed is determined by the patient's anatomy and the surgeon's preference. (see step 8).

Using only one hand to compress, the surgeon can tell when the implant is properly positioned under fluoro by watching the self-retaining washer tilt and lay flat against the bone, indicating that the device is fully compressed. The amount of compression needed is verified using “continuous fluoroscopy” and tactile feedback.

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